

Accessory: Treatment programs

Mode	Body Part	Program	Pulse rate (Hz)	Pulse width (uS)	Treatment time (Min)	Type of waveform
TENS	NECK	P1	80-120	120-100	Default:30 Adjustable:(5-90)	Modulation
		P2	4	150-200	Default:30 Adjustable:(5-90)	Modulation
		U1	Default:35 Adjustable:(2-100)	Default:200 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Continue
	SHOULDER	P1	80-100	100	Default:30 Adjustable:(5-90)	Modulation
		P2	2-60	260-160	Default:30 Adjustable:(5-90)	Modulation
		P3	Default:100 Adjustable:(2-100)	Default:150 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Burst
	ARM	P1	2	250	Default:30 Adjustable:(5-90)	Continue
		P2	100	150	Default:30 Adjustable:(5-90)	Burst
		U1	Default:100 Adjustable:(2-100)	Default:200 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Continue
	HAND	P1	100	100	Default:30 Adjustable:(5-90)	Continue
		P2	2-10	200	Default:30 Adjustable:(5-90)	Modulation
		P3	Default:60 Adjustable:(2-100)	Default:260 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Modulation
	BACK	P1	60/50/45/10/50/35	200	Default:30 Adjustable:(5-90)	Modulation
		P2	6/8/10	250	Default:30 Adjustable:(5-90)	Modulation
		U1	Default:50 Adjustable:(2-100)	Default:100 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Acupuncture
	ABDOMEN	P1	80-120	120-100	Default:30 Adjustable:(5-90)	Modulation
		P2	120	55	Default:30 Adjustable:(5-90)	Continue
		U1	Default:80 Adjustable:(2-100)	Default:100 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Continue
	HIP	P1	100	150	Default:30 Adjustable:(5-90)	Burst
		P2	40/6/50	200	Default:30 Adjustable:(5-90)	Modulation
		U1	Default:80 Adjustable:(2-100)	Default:180 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Continue
	LEG	P1	40/6/50	250	Default:30 Adjustable:(5-90)	Modulation
		P2	80	150	Default:30 Adjustable:(5-90)	Modulation
		U1	Default:6-10 Adjustable:(2-100)	Default:200 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Modulation
	FOOT	P1	80-120	100-120	Default:30 Adjustable:(5-90)	Modulation
		P2	2-10	200	Default:30 Adjustable:(5-90)	Modulation
		U1	Default:2-60 Adjustable:(2-100)	Default:260-160 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Modulation
	JOINT	P1	100	150	Default:30 Adjustable:(5-90)	Burst
		P2	120	100-120	Default:30 Adjustable:(5-90)	Modulation
		U1	Default:80 Adjustable:(2-100)	Default:180 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Continue

Accessory: Treatment programs

Mode	Body Part	Program	Pulse rate (Hz)	Pulse width (uS)	Treatment time (Min)	Type of waveform
EMS	NECK	P1	30	200	Default:30 Adjustable:(5-90)	Synchronous
		P2	40	200	Default:30 Adjustable:(5-90)	Synchronous
		U1	Default:50 Adjustable:(20-100)	Default:200 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Synchronous
	SHOULDER	P1	45	200	Default:30 Adjustable:(5-90)	Synchronous
		P2	55	200	Default:30 Adjustable:(5-90)	Synchronous
		U1	Default:80 Adjustable:(20-100)	Default:200 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Synchronous
	ARM	P1	50	150	Default:30 Adjustable:(5-90)	Synchronous
		P2	60	150	Default:30 Adjustable:(5-90)	Synchronous
		U1	Default:80 Adjustable:(20-100)	Default:150 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Synchronous
	Hand	P1	4	200	Default:30 Adjustable:(5-90)	Synchronous
		P2	5	300	Default:30 Adjustable:(5-90)	Synchronous
		U1	Default:20 Adjustable:(20-100)	Default:150 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Synchronous
	BACK	P1	60	200	Default:30 Adjustable:(5-90)	Synchronous
		P2	70	200	Default:30 Adjustable:(5-90)	Synchronous
		U1	Default:80 Adjustable:(20-100)	Default:200 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Synchronous
	ABDOMEN	P1	20	200	Default:30 Adjustable:(5-90)	Synchronous
		P2	50	200	Default:30 Adjustable:(5-90)	Synchronous
		U1	Default:60 Adjustable:(20-100)	Default:200 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Synchronous
	HIP	P1	30	150	Default:30 Adjustable:(5-90)	Synchronous
		P2	60	150	Default:30 Adjustable:(5-90)	Synchronous
		U1	Default:40 Adjustable:(20-100)	Default:150 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Synchronous
	LEG	P1	20	200	Default:30 Adjustable:(5-90)	Synchronous
		P2	80	200	Default:30 Adjustable:(5-90)	Synchronous
		U1	Default:25 Adjustable:(20-100)	Default:200 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Synchronous
	FOOT	P1	4	200	Default:30 Adjustable:(5-90)	Synchronous
		P2	5	300	Default:30 Adjustable:(5-90)	Synchronous
		U1	Default:20 Adjustable:(20-100)	Default:200 Adjustable:(100-300)	Default:30 Adjustable:(5-90)	Synchronous
MASSAGE	KNEAD	P1	28-44	120~250	30	Modulation
	RUB	P1	25-79	120~250	30	Modulation
	TAP	P1	49-97	100~240	30	Modulation

Manufactured for Caremax Australia by:

Le Reve Healthcare Pty Ltd AKA Caremax

Factory 28, 25-39 Cook Road

Mitcham Victoria Australia

Web: www.caremax.com.au

Email: enquiries@caremax.com.au